



Grace Bio-Labs  
Laboratory Method

SYPRO<sup>®</sup> Ruby Staining of ONCYTE<sup>®</sup> Porous  
Nitrocellulose Film Slides

**1.0 Purpose:**

- 1.1 General staining of proteins deposited on porous nitrocellulose films for sample protein quantification.

**2.0 Equipment & Materials**

2.1 Equipment

- 2.1.1 Orbital shaker
- 2.1.2 Microarray Fluorescence Scanner such as the GenePix 4000B Microarray Scanner (Molecular Devices)

2.2 Materials

- 2.2.1 Glass staining dishes with lids (Wheaton Cat# 900203).
- 2.2.2 Glass slide racks with handles (Wheaton Cat# 900204, 900205).
- 2.2.3 ProPlate<sup>™</sup> modules (Grace Bio-Labs)

2.3 Reagents

- 2.3.1 SYPRO<sup>®</sup> Ruby (Invitrogen Cat# S12000)
- 2.3.2 Super G Blocking Reagent, optional (Grace Bio-Labs)
- 2.3.3 Fixative Solution: 7% acetic acid, 10% methanol in distilled water
- 2.3.4 Distilled water

**3.0 Procedure**



- 3.1 Optional: Submerge spotted protein array in Super G blocking reagent for 1 hour with no agitation.
  - 3.1.1 Use of Super G blocking reagent will provide lower background fluorescence levels.
- 3.2 Transfer slides to a staining jar containing fresh water.
  - 3.2.1 Wash 3 times for 15 minutes each, with agitation (105 rpm on orbital shaker).
- 3.3 Transfer slides to a staining jar containing the fixative solution.

- 3.3.1 Incubate for 15 minutes, with agitation (105 rpm on orbital shaker).
- 3.4 Transfer slides to a staining jar containing fresh water.
  - 3.4.1 Wash 4 times for 5 minutes each, with agitation (105 rpm on orbital shaker).
- 3.5 Transfer slides to a staining jar containing the SYPRO<sup>®</sup> Ruby stain.
  - 3.5.1 Alternatively: To save on reagent, instead of a staining jar use a ProPlate<sup>™</sup> module such as the 1-well module.
  - 3.5.2 Apply a 1-well ProPlate<sup>™</sup> module to the wet slide (work quickly to prevent slide drying) and add 2 ml of SYPRO<sup>®</sup> Ruby stain.
- 3.6 Perform staining for 30 minutes at room temperature, with agitation (105 rpm on orbital shaker).
  - 3.6.1 If using ProPlate<sup>™</sup> module, perform incubation with agitation at 50 rpm.
- 3.7 Transfer slides to a staining jar containing fresh water.
  - 3.7.1 Wash 4 times for 1 minute each, with agitation (105 rpm on orbital shaker).
  - 3.7.2 If using a ProPlate<sup>™</sup> module, remove SYPRO<sup>®</sup> Ruby from chamber with a pipette and replace with 2 ml fresh water. Perform the first water wash in the ProPlate<sup>™</sup> module, then unclip and transfer slide to a staining jar for washes 2-4.
- 3.8 Air-dry the slides in the dark.
- 3.9 Scan slides using a fluorescent scanner (these settings serve as a starting point and may vary from scanner to scanner):
  - 3.9.1 Axon 4000B Scanner
    - 3.9.1.1 532 nm: 33% laser power, 400 PMT setting
    - 3.9.1.2 635 nm: 100% laser power, 800 PMT setting
- 3.10 See Appendix for typical results obtained with this method.

#### 4.0 References:

- 4.3 SYPRO<sup>®</sup> Ruby Protein Gel Stain Product Information, Molecular Probes, Rev. 11/20/2007.
- 4.4 SYPRO<sup>®</sup> Ruby Protein Stains Instruction Manual, BioRad, 4006173 Rev B.

## Appendix

**Figure 1.** SYPRO<sup>®</sup> Ruby Protein Quantification on ONCYTE<sup>®</sup> Film Slides.**Figure 2.** Fluorescent Signal from SYPRO<sup>®</sup> Ruby Stained Protein on ONCYTE<sup>®</sup> Film Slides.**Figure 2.** Fluorescent signal after staining with SYPRO<sup>®</sup> Ruby is linear over 3 orders of magnitude. Shown is the dynamic range of protein staining with SYPRO<sup>®</sup> Ruby on a SuperNOVA nitrocellulose film slide blocked with Super G blocking reagent. BSA standard was spotted in various RPPA lysate buffers with  $R^2$  averaging 0.9955 for all of the buffers (Buffer A: SDS/Tris, Buffer B: SDS/Tris/Triton; Buffer C: Urea; Buffer D: Grace Bio-Labs RPPA lysis buffer).