


SecureSlip™ Coverglass for Cell Culture



Coverglass is No. 1.5 Desag glass.

Each coverglass is supported by clear CultureWell silicone which is 0.5mm thick.

Custom sizes are available - please contact us at 1-800-813-7339.



Each pouch contains 15- SecureSlips (as shown)
Each package contains 5 pouches (5X15) = 75

MS-12 MS-12

STERILE

The diagram shows a 3x5 grid of square coverglass units. Dimensions are indicated: 17 mm for the width of the grid and 12 mm for the height of the grid. A 12 mm dimension is also shown for the width of individual units. Labels include: 'No. 1.5 coverglass (desag)', '0.50 mm clear silicone', 'Silicone Coverglass', and 'Clear plastic sheet'.

Item No.	Coverglass Size (mm)	Quantity Per Pack	Price Per Pack (US)
MS-12	12 x 12	75	\$ 123.00

Contents of each pouch are sterile.

MSR-12 MSR-12

STERILE

The diagram shows a 3x5 grid of circular coverglass units. Dimensions are indicated: 20 mm for the height of the grid and 12 mm for the height of individual units. Labels include: 'No. 1.5 coverglass (desag)', '0.50 mm clear silicone', 'Silicone Coverglass', and 'Clear plastic sheet'.

Item No.	Coverglass Size (mm)	Quantity Per Pack	Price Per Pack (US)
MSR-12	12 dia.	75	\$ 135.00

Contents of each pouch are sterile.