

# SecureSeal™ Imaging Spacers




Thin, double sided adhesive spacers confine specimens without compression.

Don't see what you need listed? Send a measured drawing and specifications to us. We will assist you in fabricating the ideal spacer material and thickness for your application


SecureSeal adhesive is now available in sheets. Sheets may be purchased in a single layer thickness or prepared in multiple layers for your convenience.

If you would like us to laser cut this material for you - please contact tech support.

SS1X9		Spacer Thickness (mm)	Quantity Per Pack	Price Per Pack (US)
 <p>— — ← 0.15 mm / without liners</p>		0.12	100	\$ 101.00

SS1X13	Item No.	Spacer Thickness (mm)	Quantity Per Pack	Price Per Pack (US)
 <p>— — ← 0.15 mm / without liners</p>	JTR13R-0.5	0.12	100	\$110.00

SS1X20		Spacer Thickness (mm)	Quantity Per Pack	Price Per Pack (US)
 <p>— — ← 0.15 mm / without liners</p>		0.12	100	\$ 124.00

SS8X9		Spacer Thickness (mm)	Quantity Per Pack	Price Per Pack (US)
 <p>— — — — — ← 0.15 mm / without liners</p>		0.12	100	\$ 141.00

## SecureSeal™ Adhesive Sheets

The same SecureSeal adhesive used to make HybriWell and SecureSeal incubation chambers. Thin, double-sided tape on a non-woven carrier.

Clear, polyester liner on one surface, white clay coated paper on the other.

Sheet size = 7.25" x 10"

1 Sheet per Pack

SA-S-	Item No.	Material Thickness (mm)	Quantity Per Pack	Price Per Pack (US)
	SA-S-1L	0.12	1	\$ 25.00
	SA-S-2L	0.24	1	\$ 36.00
	SA-S-3L	0.36	1	\$ 52.00
	SA-S-4L	0.48	1	\$ 103.00