

## SecureSeal™ Adhesive Sheets



SecureSeal adhesive is now available in sheets. Sheets may be purchased in a single layer thickness or prepared in multiple layers for your convenience.

Thin, double sided adhesive sheets allow users to personally cut spacers or gaskets to size.

If you would like us to laser cut this material for you - please contact tech support.

---

### SecureSeal™ Adhesive Sheets

SA-S-

The same SecureSeal adhesive used to make HybriWell and SecureSeal incubation chambers. Thin, double-sided tape on a non-woven carrier.

Clear, polyester liner on one surface, white clay coated paper on the other.

Sheet size = 7.25" x 10"

1 Sheet per Pack

| Item No. | Material Thickness (mm) | Quantity Per Pack | Price Per Pack (US) |
|----------|-------------------------|-------------------|---------------------|
| SA-S-1L  | 0.12                    | 1                 | \$ 25.00            |
| SA-S-2L  | 0.24                    | 1                 | \$ 36.00            |
| SA-S-3L  | 0.36                    | 1                 | \$ 52.00            |
| SA-S-4L  | 0.48                    | 1                 | \$ 103.00           |

---