

CoverWell™ Perfusion Chambers



Form removable, water-tight multiwell chambers by pressing to coverglass or microscope slides.
Thin (0.18 mm / 0.007") polycarbonate cover has excellent transparency.

PC1R							PC1R-
		Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)	
<p>25 20 25 25 1.6mm diameter ports, 2 places</p>							
<p>0.18 mm polycarbonate 0.5 mm silicone</p>	PC1R-0.5	0.6	180	40	\$ 98.00		
<p>0.18 mm polycarbonate 1.0 mm silicone</p>	PC1R-1.0	0.9	300	40	\$ 104.00		
<p>0.18 mm polycarbonate 2.0 mm silicone</p>	PC1R-2.0	1.7	520	20	\$ 65.00		
<p>0.18 mm polycarbonate 2.5 mm silicone</p>	PC1R-2.5	2.5	750	20	\$ 70.00		
PC1L							PC1L-
		Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)	
<p>43 32 3 26 19 Non adhesive peel tab 1.6mm diameter ports, 2 places</p>							
<p>0.18 mm polycarbonate 0.5 mm silicone</p>	PC1L-0.5	0.6	350	40	\$ 146.00		
<p>0.18 mm polycarbonate 1.0 mm silicone</p>	PC1L-1.0	0.9	550	40	\$ 174.00		
<p>0.18 mm polycarbonate 2.0 mm silicone</p>	PC1L-2.0	1.7	1000	20	\$ 94.00		
<p>0.18 mm polycarbonate 2.5 mm silicone</p>	PC1L-2.5	2.5	1500	20	\$ 94.00		
PC3L							PC3L-
		Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)	
<p>43 32 3 26 6.35 3 Non adhesive peel tab 1.6mm diameter ports, 6 places</p>							
<p>0.18 mm polycarbonate 0.5 mm silicone</p>	PC3L-0.5	0.6	75	40	\$ 101.00		
<p>0.18 mm polycarbonate 1.0 mm silicone</p>	PC3L-1.0	0.9	120	40	\$ 175.00		
<p>0.18 mm polycarbonate 2.0 mm silicone</p>	PC3L-2.0	1.7	200	20	\$ 96.00		
<p>0.18 mm polycarbonate 2.5 mm silicone</p>	PC3L-2.5	2.5	300	20	\$ 99.00		

PC4L		Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
PC4L-0.5		PC4L-0.5	0.6	70	40	\$158.00
PC4L-1.0		PC4L-1.0	0.9	110	40	\$187.00
PC4L-2.0		PC4L-2.0	1.7	200	20	\$ 101.00
PC4L-2.5		PC4L-2.5	2.5	300	20	\$ 107.00

PC8R		Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
PC8R-0.5		PC8R-0.5	0.6	35	20	\$ 96.00
PC8R-1.0		PC8R-1.0	0.9	60	20	\$ 108.00
PC8R-2.0		PC8R-2.0	1.7	100	10	\$ 65.00
PC8R-2.5		PC8R-2.5	2.5	150	10	\$ 70.00