

SecureSeal™ Hybridization Chambers



Form removable, peel-and-stick enclosures to isolate specimens affixed to glass microscope slides.

Each package includes illustrated instructions for use and 120 adhesive port seals. Additional seal tabs are sold separately - (listed below).

SecureSeals are recommended for most protein and nucleic acid assays. If you are using Cy5® or Alexa Fluor® 647 direct-labeled DNA probes please refer to our new “fluorescent friendly” chambers or contact tech support 1-800-813-7339 for details.

		SA50	SA50	SA50
		Item No.	Quantity Per Pack	Price Per Pack (US)
SA50		SA50	50	\$ 166.00
SA20		SA20-0.5	40	\$ 117.00
		SA20-1.0	40	\$ 125.00
		SA20-2.0	40	\$ 145.00
		SA20-2.5	40	\$ 156.00
SA8R		SA8R-0.5	20	\$ 114.00
		SA8R-1.0	20	\$ 124.00
		SA8R-2.0	20	\$ 156.00
		SA8R-2.5	20	\$ 158.00

More >>

<p>SA200</p> <p>1.6mm diameter ports, 2 places</p> <p>0.25mm polycarbonate 0.5 mm silicone adhesive</p>	<p>SA200</p> <table border="0"> <thead> <tr> <th style="text-align: left;">Item No.</th> <th style="text-align: left;">Quantity Per Pack</th> <th style="text-align: left;">Price Per Pack (US)</th> </tr> </thead> <tbody> <tr> <td>SA200</td> <td>50</td> <td>\$ 191.00</td> </tr> </tbody> </table>	Item No.	Quantity Per Pack	Price Per Pack (US)	SA200	50	\$ 191.00
Item No.	Quantity Per Pack	Price Per Pack (US)					
SA200	50	\$ 191.00					

<p>SA500</p> <p>1.6mm diameter ports, 2 places</p> <p>0.25mm polycarbonate 0.5 mm silicone adhesive</p>	<p>SA500</p> <table border="0"> <thead> <tr> <th style="text-align: left;">Item No.</th> <th style="text-align: left;">Quantity Per Pack</th> <th style="text-align: left;">Price Per Pack (US)</th> </tr> </thead> <tbody> <tr> <td>SA500</td> <td>50</td> <td>\$ 196.00</td> </tr> </tbody> </table>	Item No.	Quantity Per Pack	Price Per Pack (US)	SA500	50	\$ 196.00
Item No.	Quantity Per Pack	Price Per Pack (US)					
SA500	50	\$ 196.00					

<p>SA2657</p> <p>1.6mm diameter ports, 2 places</p> <p>0.25mm polycarbonate 0.5 mm silicone adhesive</p>	<p>SA2657</p> <table border="0"> <thead> <tr> <th style="text-align: left;">Item No.</th> <th style="text-align: left;">Quantity Per Pack</th> <th style="text-align: left;">Price Per Pack (US)</th> </tr> </thead> <tbody> <tr> <td>SA2657</td> <td>25</td> <td>\$ 184.00</td> </tr> </tbody> </table>	Item No.	Quantity Per Pack	Price Per Pack (US)	SA2657	25	\$ 184.00
Item No.	Quantity Per Pack	Price Per Pack (US)					
SA2657	25	\$ 184.00					

<p>SA2260</p> <p>1.6mm diameter ports, 2 places</p> <p>0.25mm polycarbonate 0.5 mm silicone adhesive</p>	<p>SA2260</p> <table border="0"> <thead> <tr> <th style="text-align: left;">Item No.</th> <th style="text-align: left;">Quantity Per Pack</th> <th style="text-align: left;">Price Per Pack (US)</th> </tr> </thead> <tbody> <tr> <td>SA2260</td> <td>25</td> <td>\$ 189.00</td> </tr> </tbody> </table>	Item No.	Quantity Per Pack	Price Per Pack (US)	SA2260	25	\$ 189.00
Item No.	Quantity Per Pack	Price Per Pack (US)					
SA2260	25	\$ 189.00					

<p>SA2171</p> <p>1.6mm diameter ports, 2 places</p> <p>0.25mm polycarbonate 0.5 mm silicone adhesive</p>	<p>SA2171</p> <table border="0"> <thead> <tr> <th style="text-align: left;">Item No.</th> <th style="text-align: left;">Quantity Per Pack</th> <th style="text-align: left;">Price Per Pack (US)</th> </tr> </thead> <tbody> <tr> <td>SA2171</td> <td>25</td> <td>\$ 192.00</td> </tr> </tbody> </table>	Item No.	Quantity Per Pack	Price Per Pack (US)	SA2171	25	\$ 192.00
Item No.	Quantity Per Pack	Price Per Pack (US)					
SA2171	25	\$ 192.00					

<p>SA4545</p> <p>1.6mm diameter ports, 2 places</p> <p>0.25mm polycarbonate 0.5 mm silicone adhesive</p>	<p>SA4545</p> <table border="0"> <thead> <tr> <th style="text-align: left;">Item No.</th> <th style="text-align: left;">Quantity Per Pack</th> <th style="text-align: left;">Price Per Pack (US)</th> </tr> </thead> <tbody> <tr> <td>SA4545</td> <td>20</td> <td>\$ 140.00</td> </tr> </tbody> </table>	Item No.	Quantity Per Pack	Price Per Pack (US)	SA4545	20	\$ 140.00
Item No.	Quantity Per Pack	Price Per Pack (US)					
SA4545	20	\$ 140.00					

<p>SA2X2121</p> <p>1.6mm diameter ports, 2 places</p> <p>0.25mm polycarbonate 0.5 mm silicone adhesive</p>	<p>SA2X2121</p> <table border="0"> <thead> <tr> <th style="text-align: left;">Item No.</th> <th style="text-align: left;">Quantity Per Pack</th> <th style="text-align: left;">Price Per Pack (US)</th> </tr> </thead> <tbody> <tr> <td>SA2X2121</td> <td>25</td> <td>\$ 189.00</td> </tr> </tbody> </table> <p>More >></p>	Item No.	Quantity Per Pack	Price Per Pack (US)	SA2X2121	25	\$ 189.00
Item No.	Quantity Per Pack	Price Per Pack (US)					
SA2X2121	25	\$ 189.00					

SA2X2134

75.5
34.5
2
25.5
21.5

1.6mm diameter ports, 2 places

0.25mm polycarbonate
0.5 mm silicone adhesive

SA2X2134

Item No.	Quantity Per Pack	Price Per Pack (US)
SA2X2134	25	\$ 195.00

package includes 120 adhesive port seals and illustrated instructions

SA3S-0.5

75
21
23.9
25
20

1.6mm diameter ports, 6 places

0.25mm polycarbonate
0.5 mm silicone adhesive

SA3S-0.5

Item No.	Quantity Per Pack	Price Per Pack (US)
SA3S-0.5	25	\$ 200.00

SA6L-0.5

75
9.8
11.8
25
20

1.6mm diameter ports, 12 places

0.25mm polycarbonate
0.5 mm silicone adhesive

SA6L-0.5

Item No.	Quantity Per Pack	Price Per Pack (US)
SA6L-0.5	25	\$ 205.00

SA8S

43
7
9
26
7

1.6mm diameter ports, 16 places

SA8S-

Item No.	Quantity Per Pack	Price Per Pack (US)
SA8S-0.5	25	\$ 114.00
SA8S-1.0	25	\$ 119.00

SA8S-0.5

0.25mm polycarbonate
0.5 mm silicone adhesive

SA8S-1.0

0.25mm polycarbonate
1.0 mm silicone adhesive

SA16S-0.5

75
7
9
25
7

1.6mm diameter ports, 32 places

0.25mm polycarbonate
0.5 mm silicone adhesive

SA16S-0.5

Item No.	Quantity Per Pack	Price Per Pack (US)
SA16S-0.5	25	\$ 222.00

Seal Tabs

3.8
7.62 dia
perforation
3.8
49.5
57
114

Adhesive Seal Tabs

Item No.	Quantity Per Pack	Price Per Pack (US)
ST200	200	\$ 47.00

cover filling ports to prevent evaporation