



EXPL6

Silicone Gasket Dimensions (mm)




Silicone Depth = 0.25mm

Cover Glass Dimensions



Coverglass Thickness = No. 1.5

Microscope Slide Dimensions



Slide thickness = 1mm thickness

Fisher Scientific - Tissue Path SuperFrost Plus
Part # 22-035813

ExpressLane EXPL6

Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
0.25	18	25	\$ 210.00