

HybriWell seals securely to a microscope slide to eliminate evaporation when using a small reagent volume.

Each package includes illustrated instructions for use, applicator and 200 adhesive port seals. Additional seal tabs are sold separately (listed below).

Need more or less volume? Each of the designs listed may be easily customized to increase or decrease chamber depth.


HybriWells are recommended for most protein and nucleic acid assays. If you are using Cy5® or Alexa Fluor® 647 direct-labeled DNA probes please refer to our new “fluorescent friendly” chambers (listed below) or contact tech support 1-800-813-7339 for details.


F.Y.I - All the current designs may be fabricated using our “fluorescent friendly” adhesive material.


Don't see what you need listed? Either, refer to the custom HybriWell chamber questionnaire, or send us a sketch and specifications. We will fabricate the ideal chamber for your application.

Model	Diagram	Notes	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
HBW13		Non adhesive peel tab 1.6mm diameter ports, 2 places	0.15	18	100	\$ 196.00


Model	Diagram	Notes	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
HBW20		Non adhesive peel tab 1.6mm diameter ports, 2 places	0.15	30	100	\$ 193.00

Model	Diagram	Notes	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
HBW1932		Non adhesive peel tab 1.6mm diameter ports, 2 places	0.15	30-50	100	\$ 196.00

Model	Diagram	Notes	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
HBW2222		Non adhesive peel tab 1.6mm diameter ports, 2 places	0.15	30-50	100	\$ 191.00

Model	Diagram	Notes	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
HBW75		Non adhesive peel tab 3.2mm diameter ports, 2 places	0.15	50-100	100	\$ 196.00


HBW2240



HBW2240

Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
0.25	180-200	100	\$ 202.00


HBW2260



HBW2260

Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
0.25	280-330	100	\$ 248.00


HBW2265



HBW2265

Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
0.25	310-425	100	\$ 276.00


HBW4545



HBW4545

Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
0.15	150-300	50	\$ 262.00


HBW2160



HBW2160

Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
0.25	275-320	100	\$ 228.00


HBW6021



HBW6021

Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
0.25	275-320	100	\$ 228.00

HBW3S



HBW3S

Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
0.25	100	100	\$ 257.00

<p>HBW3S</p>	<p>HBW3S</p> <table border="0"> <tr> <td>Chamber Depth (mm)</td> <td>Volume Per Well (ul)</td> <td>Quantity Per Pack</td> <td>Price Per Pack (US)</td> </tr> <tr> <td>0.25</td> <td>100</td> <td>100</td> <td>\$ 257.00</td> </tr> </table>	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)	0.25	100	100	\$ 257.00
Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)						
0.25	100	100	\$ 257.00						

<p>HBW6L</p>	<p>HBW6L</p> <table border="0"> <tr> <td>Chamber Depth (mm)</td> <td>Volume Per Well (ul)</td> <td>Quantity Per Pack</td> <td>Price Per Pack (US)</td> </tr> <tr> <td>0.25</td> <td>50</td> <td>100</td> <td>\$ 255.00</td> </tr> </table>	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)	0.25	50	100	\$ 255.00
Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)						
0.25	50	100	\$ 255.00						

<p>Seal Tabs</p>	<p>Adhesive Seal Tabs</p> <table border="0"> <tr> <td>Item No.</td> <td>Quantity Per Pack</td> <td>Price Per Pack (US)</td> </tr> <tr> <td>ST200</td> <td>200</td> <td>\$ 47.00</td> </tr> </table> <p>cover filling ports to prevent evaporation</p>	Item No.	Quantity Per Pack	Price Per Pack (US)	ST200	200	\$ 47.00
Item No.	Quantity Per Pack	Price Per Pack (US)					
ST200	200	\$ 47.00					

Fluorescent Friendly Chambers for Sensitive DNA Microarrays



These new microarray chambers are designed specifically for compatibility with Cy5[®] or AlexaFluor[®] 647 direct labeled DNA probes for hybridization assays. These chambers as well as the original HybriWells may also be used with all other fluorochrome labeled probes.

Each package includes illustrated instructions for use, applicator and 200 adhesive port seals. Additional seal tabs are sold separately (listed below).


Need more or less volume? Each of the designs listed may be easily customized to increase or decrease chamber depth.

F.Y.I - All the current HybriWell designs listed above may be fabricated using our “fluorescent friendly” adhesive material.

Don't see what you need listed? Either, refer to the custom HybriWell chamber questionnaire, or send us a sketch and specifications. We will fabricate the ideal chamber for your application.

<p>HBW2222-FL</p>	<p>HBW2222-FL</p> <table border="0"> <tr> <td>Chamber Depth (mm)</td> <td>Volume Per Well (ul)</td> <td>Quantity Per Pack</td> <td>Price Per Pack (US)</td> </tr> <tr> <td>0.25</td> <td>30-50</td> <td>100</td> <td>\$ 221.00</td> </tr> </table>	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)	0.25	30-50	100	\$ 221.00
Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)						
0.25	30-50	100	\$ 221.00						


HBW2240-FL



HBW2240-FL

Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
0.25	180-200	100	\$ 226.00

HBW2260-FL



HBW2260-FL

Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
0.25	280-330	100	\$ 245.00