

# HybriSlip™ Hybridization Covers



HybriSlips are manufactured using 0.25mm polycarbonate for rigidity and temperature resistance.

HybriSlips are packaged in quantities of 100 per bag in strips of 5 each, with protective liner on both surfaces.

Popular sizes are listed below. If you don't see what you need listed, HybriSlips are easily cut to custom sizes.

HS22		HS22 / HS22-CS									
	<table border="0"> <thead> <tr> <th style="text-align: left;">Item No.</th> <th style="text-align: left;">Quantity Per Pack</th> <th style="text-align: left;">Price Per Pack (US)</th> </tr> </thead> <tbody> <tr> <td>HS22</td> <td>100</td> <td>\$ 22.00</td> </tr> <tr> <td>HS22-CS</td> <td>1000</td> <td>\$ 191.00</td> </tr> </tbody> </table>	Item No.	Quantity Per Pack	Price Per Pack (US)	HS22	100	\$ 22.00	HS22-CS	1000	\$ 191.00	
Item No.	Quantity Per Pack	Price Per Pack (US)									
HS22	100	\$ 22.00									
HS22-CS	1000	\$ 191.00									


Hs2525		HS2525 / HS2525-CS									
	<table border="0"> <thead> <tr> <th style="text-align: left;">Item No.</th> <th style="text-align: left;">Quantity Per Pack</th> <th style="text-align: left;">Price Per Pack (US)</th> </tr> </thead> <tbody> <tr> <td>Hs2525</td> <td>100</td> <td>\$ 30.00</td> </tr> <tr> <td>HS2525-CS</td> <td>1000</td> <td>\$ 269.00</td> </tr> </tbody> </table>	Item No.	Quantity Per Pack	Price Per Pack (US)	Hs2525	100	\$ 30.00	HS2525-CS	1000	\$ 269.00	
Item No.	Quantity Per Pack	Price Per Pack (US)									
Hs2525	100	\$ 30.00									
HS2525-CS	1000	\$ 269.00									

HS40		HS40 / HS40-CS									
	<table border="0"> <thead> <tr> <th style="text-align: left;">Item No.</th> <th style="text-align: left;">Quantity Per Pack</th> <th style="text-align: left;">Price Per Pack (US)</th> </tr> </thead> <tbody> <tr> <td>HS40</td> <td>100</td> <td>\$ 26.00</td> </tr> <tr> <td>HS40-CS</td> <td>1000</td> <td>\$ 235.00</td> </tr> </tbody> </table>	Item No.	Quantity Per Pack	Price Per Pack (US)	HS40	100	\$ 26.00	HS40-CS	1000	\$ 235.00	
Item No.	Quantity Per Pack	Price Per Pack (US)									
HS40	100	\$ 26.00									
HS40-CS	1000	\$ 235.00									

HS4024		HS4024 / HS4024-CS									
	<table border="0"> <thead> <tr> <th style="text-align: left;">Item No.</th> <th style="text-align: left;">Quantity Per Pack</th> <th style="text-align: left;">Price Per Pack (US)</th> </tr> </thead> <tbody> <tr> <td>HS4024</td> <td>100</td> <td>\$ 30.00</td> </tr> <tr> <td>HS4024-CS</td> <td>1000</td> <td>\$ 268.00</td> </tr> </tbody> </table>	Item No.	Quantity Per Pack	Price Per Pack (US)	HS4024	100	\$ 30.00	HS4024-CS	1000	\$ 268.00	
Item No.	Quantity Per Pack	Price Per Pack (US)									
HS4024	100	\$ 30.00									
HS4024-CS	1000	\$ 268.00									

HS60

HS60 / HS60-CS




← 0.25 mm polycarbonate

Item No.	Quantity Per Pack	Price Per Pack (US)
HS60	100	\$ 30.00
HS60-CS	1000	\$ 269.00

HS6024

HS6024 / HS6024-CS




← 0.25 mm polycarbonate

Item No.	Quantity Per Pack	Price Per Pack (US)
HS6024	100	\$ 35.00
HS6024-CS	1000	\$ 308.00

HS4550

HS4550 / HS4550-CS



← 0.25 mm polycarbonate

Item No.	Quantity Per Pack	Price Per Pack (US)
HS4550	50	\$ 54.00
HS4550-CS	250	\$ 480.00