


ProPlate™ Multi-Array Slide System



ProPlate Slide Module



204862


Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
15 - 350	2	\$ 49.00

Stainless Steel Spring Clips

204838 / 204839



Spring Clip End View



Slide module assembled using Spring Clips End View

Slide module assembled using Spring Clips Top View


Stainless steel spring clips are recommended for use with glass slides which have un-polished edges.

Item No.	Quantity Per Pack	Price Per Pack (US)
204838 (un-numbered)	2	\$ 22.00
204839 * (with numbering)	2	\$ 27.00

* Numbered clips are permanently laser etched along the top surface, marking wells at #1-15 left clip (top), # 2-16 on the right clip (bottom). Custom printing is also available.

Seal Strips and Applicator


204874



Item No.	Quantity Applicators Per Pack	Quantity Seal Strips Per Pack	Price Per Pack (US)
204874	1	50	\$ 35.00

ProPlate Tray Set

204860



Each ProPlate Tray Set comes with;
 4, 16-well slide modules
 1 Tray and Cover
 10 Seal strips and 1 seal strip applicator.

Item No.	Quantity Per Pack	Price Per Pack (US)
204860	1	\$ 119.00

Tray conforms to SBS standards for footprint (81.7mm width x 123.7mm length) and well spacing (96 well plate).