

# Protein Crystallization Covers for 96 Well Plates




Isolate individual wells using peel & stick plate covers.

Hydrophobic, UV transparent plastic seals to plates using removable pressure sensitive adhesive.

Available as a single “un-cut” 96 well cover sheet or as a cover sheet which applies as a single sheet and allows for the removal of pre-cut individual “well caps.”

45232 / 45233

## “Un-cut” Protein Crystallization Covers




| Item No. | Liner Material | Quantity Per Pack | Price Per Pack (US) |
|----------|----------------|-------------------|---------------------|
| 45232    | Paper          | 5                 | \$ 27.00            |
| 45233*   | Plastic        | 5                 | \$ 27.00            |

\* 44233 formerly known as 4423H

45236 / 45237

## “Cut” Protein Crystallization Covers



| Item No. | Liner Material | Quantity Per Pack | Price Per Pack (US) |
|----------|----------------|-------------------|---------------------|
| 45236    | Paper          | 5                 | \$ 36.25            |
| 45237*   | Plastic        | 5                 | \$ 36.25            |

\* 45237 formerly known as 4423K