

CoverWell™ Imaging Chambers




Form removable, water-tight chambers to support thick and free-floating specimens for imaging.

Thin (0.18 mm / 0.007") polycarbonate cover has excellent transparency.

For additional sealing or permanent mounting in aqueous media, select "Adhesive Imaging Chambers"


PCI



Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
PCI-0.5	0.6	180	40	\$ 79.00
PCI-1.0	0.9	300	40	\$ 82.00
PCI-2.0	1.7	520	20	\$ 47.00
PCI-2.5	2.5	750	20	\$ 50.00

Labels for diagrams: 0.18 mm polycarbonate, 0.5 mm silicone, 1.0 mm silicone, 2.0 mm silicone, 2.5 mm silicone.


PCI-A Adhesive Imaging Chambers



Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
PCI-A-0.5	0.6	180	40	\$ 88.00
PCI-A-1.0	0.9	300	40	\$ 91.00
PCI-A-2.0	1.7	520	20	\$ 52.00
PCI-A-2.5	2.5	750	20	\$ 55.00

Labels for diagrams: 0.18 mm polycarbonate, 0.5 mm silicone adhesive, 1.0 mm silicone adhesive, 0.28 mm polycarbonate, 2.0 mm silicone adhesive, 0.18 mm polycarbonate, 2.5 mm silicone adhesive.

PCI-PAF



Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
PCI-PAF	0.26	180	25	\$ 72.50

Labels for diagram: 0.18 mm polycarbonate, 0.25 mm silicone.