

CoverWell™ Imaging Chambers



Form removable, water-tight chambers to support thick and free-floating specimens for imaging.

Thin (0.18 mm / 0.007") polycarbonate cover has excellent transparency.

For additional sealing or permanent mounting in aqueous media, select "Adhesive Imaging Chambers"

PCI					PCI-	
		Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
PCI-0.5		PCI-0.5	0.6	180	40	\$ 95.00
PCI-1.0		PCI-1.0	0.9	300	40	\$ 98.00
PCI-2.0		PCI-2.0	1.7	520	20	\$ 57.00
PCI-2.5		PCI-2.5	2.5	750	20	\$ 60.00

PCI-A	Adhesive Imaging Chambers				PCI-A-	
		Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
PCI-A-0.5		PCI-A-0.5	0.6	180	40	\$ 106.00
PCI-A-1.0		PCI-A-1.0	0.9	300	40	\$ 109.00
PCI-A-2.0		PCI-A-2.0	1.7	520	20	\$ 63.00
PCI-A-2.5		PCI-A-2.5	2.5	750	20	\$ 66.00

PCI-PAF					PCI-PAF	
		Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
		PCI-PAF	0.26	180	25	\$ 98.00