

EliSpot Slide



Convenient microscope slide format for measuring immune cell activation at the single cell level. The interior of each well is coated with a PVDF membrane.

Removable, 16 well chamber and lid is made of low protein binding material and forms a leak-free barrier between wells.


Standard slide design facilitates quick and efficient washing, handling and sample storage.

Glass surface provides optimum flatness for scanning results.


Available in glass and plastic slide types. Sold in sterile, 2-pack trays, 8 slides per order.

STERILE


EliSpot Slide



PVDF Membrane Diameter (mm)	Culture Area Per Well (cm ²)	Suggested Volume Per Well (ul)
6.0	0.4	100 - 250



Item No.	Slide Material	Quantity Per Pack	Price Per Pack (US)
42489G	Glass	8	\$ 146.00



Item No.	Slide Material	Quantity Per Pack	Price Per Pack (US)
42489P	Plastic	8	\$ 146.00