

CultureWell™ Silicone Gaskets



Isolate cells using removable hydrophobic barriers on coverglass.



Reusable CultureWell gaskets are available bulk packaged.

Gaskets are non-sterile and may be sterilized by autoclaving, UV or alcohol.

CW 2R-

NON-STERILE

CW 2R-




Item No.	Silicone Well Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
CW 2R-1.0	1.0	250 - 400	10	\$ 56.00
CW 2R-2.0	1.5	300 - 500	10	\$ 64.00

CW 3S-1.0

NON-STERILE

CW 3S-1.0




Item No.	Silicone Well Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
CW 3S-1.0	1.0	300 - 500	10	\$ 64.00

CW 4R-1.0

NON-STERILE

CW 4R-1.0




Item No.	Silicone Well Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
CW 4R-1.0	1.0	50 - 100	10	\$ 58.00

CW 8R-1.0

NON-STERILE

CW 8R-1.0




Item No.	Silicone Well Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
CW 8R-1.0	1.0	15 - 30	10	\$ 59.00

CW 50R-1.0

NON-STERILE

CW 50R-1.0



Item No.	Silicone Well Depth (mm)	Volume Per Well (ul)	Quantity Per Pack	Price Per Pack (US)
CW 50R-1.0	1.0	3 - 10	10	\$ 63.00