

# Press-to-Seal Silicone Isolators & FastWells



Isolate specimens using removable hydrophobic barriers.

Don't see what you need listed? Please refer to the Custom Silicone Isolator Questionnaire or send a measured drawing and specifications to us. We will assist you in fabricating the ideal isolator.

NEW! Each of the current seal designs is also available in CultureWell "water-clear" silicone which is highly recommend for sensitive cell culture applications.


Silicone Isolators are available in a variety of silicone thickness and where additional sealing is required, SecureSeal Adhesive on one or both surfaces is recommended.

<p>664112</p>	<p>FastWell / FW9</p>				
<table border="0"> <thead> <tr> <th style="text-align: left;">Item No.</th> <th style="text-align: left;">Quantity Per Pack</th> </tr> </thead> <tbody> <tr> <td>664112</td> <td>50</td> </tr> </tbody> </table>	Item No.	Quantity Per Pack	664112	50	
Item No.	Quantity Per Pack				
664112	50				

<p>664113</p>	<p>FastWell / FW20</p>				
<table border="0"> <thead> <tr> <th style="text-align: left;">Item No.</th> <th style="text-align: left;">Quantity Per Pack</th> </tr> </thead> <tbody> <tr> <td>664113</td> <td>50</td> </tr> </tbody> </table>	Item No.	Quantity Per Pack	664113	50	
Item No.	Quantity Per Pack				
664113	50				

<p>664114</p>	<p>FastWell / FW1932</p>				
<table border="0"> <thead> <tr> <th style="text-align: left;">Item No.</th> <th style="text-align: left;">Quantity Per Pack</th> </tr> </thead> <tbody> <tr> <td>664114</td> <td>50</td> </tr> </tbody> </table>	Item No.	Quantity Per Pack	664114	50	
Item No.	Quantity Per Pack				
664114	50				

664107-507 Non-adhesive - reusable / autoclavable



- ← 0.5 mm silicone
- ← 1.0 mm silicone
- ← 2.0 mm silicone
- ← 2.5 mm silicone


Item No.	Description	Quantity Per Pack
664507	JTR13R-0.5	25
664107	JTR13R-1.0	25
664207	JTR13R-2.0	25
664307	JTR13R-2.5	25

665107-507 Adhesive one side




665507	JTR13R-A-0.5	25
665107	JTR13R-A-1.0	25
665207	JTR13R-A-2.0	25
665307	JTR13R-A-2.5	25

666107-507 Adhesive on both sides



666507	JTR13R-A2-0.5	25
666107	JTR13R-A2-1.0	25
666207	JTR13R-A2-2.0	25
666307	JTR13R-A2-2.5	25

664104-504 Non-Adhesive - reusable / autoclavable



- ← 0.5 mm silicone
- ← 1.0 mm silicone
- ← 2.0 mm silicone
- ← 2.5 mm silicone


Item No.	Description	Quantity Per Pack	JTR20R
664504	JTR20R-0.5	25	
664104	JTR20R-1.0	25	
664204	JTR20R-2.0	25	
664304	JTR20R-2.5	25	

665104-504 SecureSeal Adhesive (PSA) one side




665504	JTR20R-A-0.5	25
665104	JTR20R-A-1.0	25
665204	JTR20R-A-2.0	25
665304	JTR20R-A-2.5	25

666104-504 SecureSeal Adhesive (PSA) both sides



666504	JTR20R-A2-0.5	25
666104	JTR20R-A2-1.0	25
666204	JTR20R-A2-2.0	25
666304	JTR20R-A2-2.5	25


664103-503 Non-adhesive - reusable / autoclavable



- ← 0.5 mm silicone
- ← 1.0 mm silicone
- ← 2.0 mm silicone
- ← 2.5 mm silicone


Item No.	Description	Quantity Per Pack	JTR1L
664503	JTR1L-0.5	25	
664103	JTR1L-1.0	25	
664203	JTR1L-2.0	25	
664303	JTR1L-2.5	25	

665103-503 Adhesive one side



665503	JTR1L-A-0.5	25
665103	JTR1L-A-1.0	25
665203	JTR1L-A-2.0	25
665303	JTR1L-A-2.5	25



666103-503 Adhesive on both sides



666503	JTR1L-A2-0.5	25
666103	JTR1L-A2-1.0	25
666203	JTR1L-A2-2.0	25
666303	JTR1L-A2-2.5	25


664101-501 Non-adhesive - reusable / autoclavable

JTR8R




Item No.	Description	Quantity Per Pack
664501	JTR8R-0.5	25
664101	JTR8R-1.0	25
664201	JTR8R-2.0	25
664301	JTR8R-2.5	25

665101-501 Adhesive one side



665501	JTR8R-A-0.5	25
665101	JTR8R-A-1.0	25
665201	JTR8R-A-2.0	25
665301	JTR8R-A-2.5	25



666105-505 Adhesive on both sides



666505	JTR8R-A2-0.5	25
666105	JTR8R-A2-1.0	25
666205	JTR8R-A2-2.0	25
666305	JTR8R-A2-2.5	25


664106-506 Non-adhesive - reusable / autoclavable

JTR12R




Item No	Description	Quantity Per Pack
664506	JTR12R-0.5	25
664106	JTR12R-1.0	25
664206	JTR12R-2.0	25

665106-506 Adhesive one side



665506	JTR12R-A-0.5	25
665106	JTR12R-A-1.0	25
665206	JTR12R-A-2.0	25


666106-506 Adhesive on both sides



666506	JTR12R-A2-0.5	25
666106	JTR12R-A2-1.0	25
666206	JTR12R-A2-2.0	25

664108-508 Non-adhesive - reusable / autoclavable


JTR24R



- ← 0.5 mm silicone
- ← 1.0 mm silicone
- ← 2.0 mm silicone
- ← 2.5 mm silicone


Item No.	Description	Quantity Per Pack
664508	JTR24R-0.5	25
664108	JTR24R-1.0	25
664208	JTR24R-2.0	25

665108-508 Adhesive one side



665508	JTR24R-A-0.5	25
665108	JTR24R-A-1.0	25
665208	JTR24R-A-2.0	25


666108-508 Adhesive on both sides



666508	JTR24R-A2-0.5	25
666108	JTR24R-A2-1.0	25
666208	JTR24R-A2-2.0	25

664115-515 Non-adhesive - reusable / autoclavable


JTR8S



- ← 0.5 mm silicone
- ← 1.0 mm silicone
- ← 2.0 mm silicone
- ← 2.5 mm silicone


Item No.	Description	Quantity Per Pack
664515	JTR8S-0.5	25
664115	JTR8S-1.0	25
664215	JTR8S-2.0	25

666115-515 Adhesive one side





665515	JTR8S-A-0.5	25
665115	JTR8S-A-1.0	25
665215	JTR8S-A-2.0	25

666115-515 Adhesive on both sides



666515	JTR8S-A2-0.5	25
666115	JTR8S-A2-1.0	25
666215	JTR8S-A2-2.0	25


664116-516 Non-adhesive - reusable / autoclavable



JTR10S


Item No.	Description	Quantity Per Pack
664516	JTR10S-0.5	25
664116	JTR10S-1.0	25
664216	JTR10S-2.0	25

665116-516 Adhesive one side





665516	JTR10S-A-0.5	25
665116	JTR10S-A-1.0	25
665216	JTR10S-A-2.0	25

666116-516 Adhesive on both sides



666516	JTR10S-A2-0.5	25
666116	JTR10S-A2-1.0	25
666216	JTR10S-A2-2.0	25


664117-517 Non-adhesive - reusable / autoclavable



JTR16S


Item No.	Description	Quantity Per Pack
664517	JTR16S-0.5	25
664117	JTR16S-1.0	25
664217	JTR16S-2.0	25

665117-517 Adhesive one side



665517	JTR16S-A-0.5	25
665117	JTR16S-A-1.0	25
665217	JTR16S-A-2.0	25

666117-517 Adhesive on both sides



666517	JTR16S-A2-0.5	25
666117	JTR16S-A2-1.0	25
666217	JTR16S-A2-2.0	25

## RED SILICONE SHEET MATERIAL

Each Sheet measures 13cm x 15cm, surfaces are protected with removable liners.

Item No.	Description	Thickness (mm)	JTR-S Sheets Per Pack
4182/283/384/581 Non-adhesive - reusable / autoclavable			
664581	JTR-S-0.5	0.5	5
664182	JTR-S-1.0	0.8	5
664283	JTR-S-2.0	1.6	5
664384	JTR-S-2.5	2.4	5
5182/283/384/581 Adhesive one side			
665581	JTR-SA-0.5	0.6	5
665182	JTR-SA-1.0	0.9	5
665283	JTR-SA-2.0	1.7	5
665384	JTR-SA-2.5	2.5	5
6182/283/384/581 Adhesive on both sides			
666581	JTR-SA2-0.5	0.8	5
666182	JTR-SA2-1.0	1.0	5
666283	JTR-SA2-2.0	1.8	5
666384	JTR-SA2-2.5	2.6	5

## CULTUREWELL SILICONE SHEET MATERIAL

Each Sheet measures 13cm x 15cm, surfaces are protected with removable liners.

Item No.	Description	Thickness (mm)	CWS-S Sheets Per Pack
Non-adhesive - reusable / autoclavable			
664475	CWS-S-0.25	0.25	5
664571	CWS-S-0.5	0.5	5
664172	CWS-S-1.0	0.8	5
664273	CWS-S-2.0	1.6	5

<< Previous