



Isolate specimens using adhesive multiwell chambers.  
 Injection molded gaskets are coated with SecureSeal PLUS adhesive layer with removable liner.

Top surface by be sealed with FlexWell Adhesive Seal Strips to prevent evaporation.

Multiple FlexWells may be inserted in the FlexWell tray for standard microtitre plate spacing to facilitate handling.

204908


FlexWell 8



Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
204908	3.3	15-150	10

204916


FlexWell 16



Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
204916	3.3	15-150	10



204064

FlexWell 64




Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
204064	1.8	3-15	5

FlexWell Seal Strips & Applicator



Item No.	Strips Per Pack	Applicators Per Pack
204974	50	1


FlexWell Slide Tray 204970



Item No.	Quantity Per Pack
204970	1

204996

FlexWell 96



Item No.	Quantity Per Pack
204996	5
204997 With Frame	1