

# HybriSlip™ Hybridization Covers




HybriSlips are manufactured using 0.25mm polycarbonate for rigidity and temperature resistance.

HybriSlips are packaged in quantities of 100 per bag in strips of 5 each, with protective liner on both surfaces.

Popular sizes are listed below. If you don't see what you need listed, HybriSlips are easily cut to custom sizes.

712222/722222

HS22 / HS22-CS




← 0.25 mm polycarbonate

Description	Item No.	Quantity Per Pack
HS22	712222	100
HS22-CS	722222	1000

712525/702525

HS25 / HS25-CS




← 0.25 mm polycarbonate

Description	Item No.	Quantity Per Pack
HS25	712525	100
HS25-CS	702525	1000

714022/724022

HS40 / HS40-CS




← 0.25 mm polycarbonate

Description	Item No.	Quantity Per Pack
HS40	714022	100
HS40-CS	724022	1000

716022/726022

HS60 / HS60-CS




← 0.25 mm polycarbonate

Description	Item No.	Quantity Per Pack
Hs60	716022	100
HS60-CS	726022	1000

716024/726024


HS2274 / HS2274-CS



Description	Item No.	Quantity Per Pack
HS2274	716024	100
HS2274-CS	726024	1000

702430/722430


HS2430 / HS2430-CS



Description	Item No.	Quantity Per Pack
HS2430	702430	100
HS2430-CS	722430	1000

714024/724024


HS4024 / HS4024-CS



Description	Item No.	Quantity Per Pack
HS4024	714024	100
HS4024-CS	724024	1000

716024/726024


HS6024 / HS6024-CS



Description	Item No.	Quantity Per Pack
HS6024	716024	100
HS6024-CS	726024	1000

714550/724550

HS4550 / HS4550-CS



Description	Item No.	Quantity Per Pack
HS4550	714550	50
HS4550-CS	724550	250