

CoverWell™ Incubation Chambers





Incubation chambers that attach without adhesive.

Reusable - easy to apply

645401/645501



PC20 / PC50



Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
645401	0.2	20	25
645501	0.5	20	50

645402/645502


PC200 / PC500



Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
645402	0.2	200	25
645502	0.5	500	50

647001/647005

PCUA / PCUA-5



Reusable holder for waterbath incubations.

Item No.	Quantity Per Pack
647001	1
647005	5