





Thin, double sided adhesive spacers confine specimens without compression.

Don't see what you need listed? Send a measured drawing and specifications to us. We will assist you in fabricating the ideal spacer material and thickness for your application

SecureSeal adhesive is now available in sheets. Sheets may be purchased in a single layer thickness or prepared in multiple layers for your convenience.

If you would like us to laser cut this material for you - please contact tech support.

Item No.	Spacer Thickness (mm)	Quantity Per Pack
654002	0.12	100
		
654004	0.12	100
		
654006	0.12	100
		
654008	0.12	100
		

## SecureSeal™ Adhesive Sheets

SA-S-

The same SecureSeal adhesive used to make HybriWell and SecureSeal incubation chambers. Thin, double-sided tape on a non-woven carrier.

Clear, polyester liner on one surface, white clay coated paper on the other.

Sheet size = 7.25" x 10"      1 Sheet per Pack

Description	Item No.	Material Thickness (mm)	Quantity Per Pack
SA-S-1L	620001	0.12	1
SA-S-2L	620002	0.24	1
SA-S-3L	620003	0.36	1
SA-S-4L	620004	0.48	1