


SecureSlip™ Coverglass for Cell Culture



Coverglass is No. 1.5 Desag glass.

Each coverglass is supported by clear CultureWell silicone which is 0.5mm thick.


Custom sizes are available - please contact us at 1-800-813-7339.



Each pouch contains 15- SecureSlips (as shown)
Each package contains 5 pouches (5X15) = 75

104112 **STERILE**

MS-12




Item No.	Coverglass Size (mm)	Quantity Per Pack
104112	12 x 12	75

Contents of each pouch are sterile.

104212 **STERILE**

MSR-12



Item No.	Coverglass Size (mm)	Quantity Per Pack
104212	12 DIA.	75

Contents of each pouch are sterile.